

# Cryopreservation of hiPSCs

---

**Equipment:**

- Cryovials

**Reagents:**

- DMSO
- iPSC Cell Culture Media (e.g. Essential 8, mTeSR, etc.)
- Freezing Media
  - Information on how to make iPSC Freezing Media is in the Preparation Step

**Procedure:****Preparation Step:**

- Label cryovials:
  - Line Name
  - Passage Number
  - Date frozen down
  - Cell Count (if applicable)
- Prepare Freezing Media (volume depends on how many vials you plan to freeze down):
  - 10% DMSO
  - 90% iPSC Cell Culture Media

**Procedure:**

1. Aspirate media from cells
2. Add 1 mL of ReleSR per well, shake to cover the well, and aspirate off the ReleSR.
3. Place plate in 37C and let incubate for 3-5 minutes
4. After cells have finished incubating, tap the plate against the hand to dislodge cells and break up cells with 1 mL of iPSC Cell Culture Media.
  - a. Be careful to not triturate cells too much. Single-cell suspensions usually have a very low recovery after thaw.
5. Spin cells down at 200g for 5 minutes
6. Aspirate the media and resuspend cells gently with 1mL of freezing media and gently break the pellet apart.
7. Add the remaining volume of Freezing Media and resuspend.
  - a. Volume to add:
    - i. Freezing ½ well per vial:
      - $(\text{\# of wells freezing down}) \times 2) - 1\text{mL} = \text{Volume to add to cells}$
    - ii. Freezing 1 well per vial
      - $(\text{\# of wells freezing down}) - 1\text{mL} = \text{Volume to add to cells}$
8. Place vials in a Mr. Frosty that's **at room temp** (not one that came right from the freezer) and place vials in a -80C freezer.
  - a. Cells will freeze at approximately 1K/min
  - b. Do not place cells in a Mr. Frosty that has come directly from a -80C freezer! Cells will freeze too quickly which will lead to cell membrane damage.

9. Leave cells in Mr. Frosty overnight and transfer cells to a Liquid Nitrogen Dewar 24 hours afterward for long-term storage.

**Protocol Updates:**

<b>Date</b>	<b>Update Notes</b>
9/19/23	Lexi – Made a change to the Cryopreservation media (now 90% E8 & 10% DMSO). Also making a change to passaging with ReleSR instead of EDTA
10/23/23	Lexi – To make the protocol more accessible for cells growing with other types of media, I changed E8 to iPSC Cell Culture media